

SAUNA BELTS

Odyssey Products is proud to provide you Sauna Belts

Sauna Belts comes in three ranges :

1. Exclusive Heating

2. Vibration & Heating

3. Magnetic, Vibration & Heating



Sauna Belt helps you sweat and relieve back pain!

Made for men and women, simply wrap the Sauna Belt around your body and with the press of a button you're sweating!

ODESSEY SAUNA BELTS

The **Sauna Belt** easily helps you to sweat, burn calories, and even eases muscle back pain! You can easily enjoy the sauna heat in the comfort of your own home!

Heat is a therapy that may in most instances allow the calming of muscles and tension in the body. Spinal complications are very hard to live with on a general everyday basis, and also make it difficult to sleep in such discomfort. Like ice packs, back heat therapy may help in muscle relaxation and even on bruises to either the skin or other tissue. The decrease of inflammation is also promoted using heat therapy.

The Sauna Belt focuses the sauna heat to help soothe back pain. The Sauna Belt is easily wrapped around your body secured with built in Velcro straps to maintain the heat on your aching back.

- Perfect for relaxation and soothing certain body aches.
- Used by professional therapists and body massagers around the world.



The Sauna Belt is a system that produces heat to help relax muscles and reduce tensions on the body. Position the Sauna Belt to the targeted area, e.g. the abdomen, waist, arms, thighs, legs, and back. Wearing the Sauna Belt each day can help relax body and soothe aches. The Amazing Brand New Sauna Belt helps you to sweat away unwanted weight, ease aches and pains in various locations, and de-stress with great convenience when and where you choose. Sauna Belt focuses the heat on the most problematic body parts to help flush out and eliminate toxins. This portable sauna is great for the abdomen, waist, back & hips. Just apply around the chosen body part and let the belt do all the work. The Sauna Belt helps you flush out toxins, reduce muscle pain, and sweat away pounds and cellulite. Using the belt is easy; simply wrap the belt around the desired area, set the temperature and enjoy. Each package comes with temperature control, diet guide, instruction manual, and a tape measure to see how many inches you have sweated off.

FEATURES

- Sweats away unwanted fat .
- You can lose weight and ease muscle pain while sitting in your own homes.
- Great for abdomen, waist, back and hips.
- Improves circulatory system of body .
- Helps maintain health and tender skin.
- Eliminates superfluous fat and beautifies figure.

PACKING DETAILS

SAUNA BELTS

Exclusive Heating

Vibration & Heating

Magnetic, Vibration & Heating

PACKING DETAILS

Quantity: 60 pieces / carton **Dimension:** 28" X 16.5" X 31" (L X W X H)

Quantity: 60 pieces / carton **Dimension:** 28" X 16.5" X 31" (L X W X H)

Quantity: 60 pieces / carton **Dimension:** 28" X 16.5" X 31" (L X W X H)

SAFETY & PRECAUTIONS

Before placing the Sauna Belt on your skin, make sure you are in a comfortable position either lying down or standing up. However, do not use the Sauna Belt if you are lying on your stomach or sitting down, as this could cause problems in the devices electrical circuit.

- Do not use the Sauna Belt longer than 50 minutes. **Remember to switch off the sauna belt after 50 minutes!**
- Put on the Sauna Belt a little loosely (**Not too tight!**), in such a way that you can check every 10 minutes if your skin is getting red. If your skin gets red, reduce the temperature to avoid possible skin injury.
- If you have **sensitive skin**, we strongly advise you to use the Sauna Belt over a T-shirt due to its heat in order to refrain from skin irritations or burns.



- It is important to use the Sauna Belt at room temperature and not to use air conditioning to try and increase its effectiveness. In this way, your body will not feel the effects of sudden changes in temperature.
- Do not use any kind of skin cream or similar product while using the Sauna Belt as your skin pores will be fully open and skin cream could cause an allergic reaction.
- Drinking hot drinks such as milk, tea or coffee while using the Sauna Belt is recommended. Do not drink cold drinks. Do not use the Sauna Belt immediately after meals; wait at least two hours after eating.

For best results, you should use the Sauna Belt for 50 minutes at a time. **Do not exceed 50 minutes per session.** It is recommended to use the sauna belt for 40 minutes. Do not use the Sauna Belt for periods of more than 45 days, wait a week before using it again.